


pH, redox and I.S.E. electrodes



Complete and clear offer which makes easy the choice of the suitable model without user's confusion.
pH electrodes: four models for portable instrument, five basic laboratory models and four models for special applications.

Every basic model is also manufactured with built-in temperature probe, T version.
The system ContATC (patented) ensures maximum response speed in the temperature measurement.

Redox and new range of last generation I.S.E. electrodes are available too.
The I.S.E. do not require internal electrolyte and their solid state membranes have long duration.



 **CRISON**
— Experts in pH —